

Prep

Make Stock Leachate

First make 400mL of each stock solution from the 1000 mg/L prepared leachate

☐ Make 10 mg/L stock:

$$1000 \text{ mg/L} * V1 = 10 \text{ mg/L} * 400\text{mL}$$

$$V1 = 4\text{mL}$$

$$400\text{mL of } 10\text{mg/L stock} = 4\text{mL of } 1000\text{mg/L stock} + 396\text{mL of FSW}$$

☐ Make 1mg/L stock:

$$10\text{mg/L} * V1 = 1\text{mg/L} * 400\text{mL}$$

$$V1 = 40\text{mL}$$

$$400\text{mL of } 1\text{mg/L stock} = 40\text{mL of } 10\text{mg/L stock} + 360\text{mL of FSW}$$

☐ Make 0.1 mg/L stock:

$$1\text{mg/L} * V1 = 0.1\text{mg/L} * 400\text{mL}$$

$$V1 = 40\text{mL}$$

$$400\text{mL of } 0.1\text{mg/L stock} = 40\text{mL of } 1\text{mg/L stock} + 360\text{mL of FSW}$$

2. Then dilute the stock to each treatment vial/jar

For a 20mL scintillation vial for coral embryo experiments with a final target volume of 15mL:

1mg/L HIGH treatment

$$10\text{mg/L} * V1 = 1\text{mg/L} * 19\text{mL}$$

$$V1 = 1.9\text{mL}$$

$$19\text{mL of } 1\text{mg/L leachate} = 17.1\text{mL of FSW} + 1.9\text{mL of } 10\text{mg/L stock}$$

0.1 mg/L MID treatment

$$1\text{mg/L} * V1 = 0.1\text{mg/L} * 19\text{mL}$$

$$V1 = 1.9\text{mL}$$

$$19\text{mL of } 0.1\text{mg/L leachate} = 17.1\text{mL of FSW} + 1.9\text{mL of } 1\text{mg/L stock}$$

0.01 mg/L LOW treatment

$$0.1\text{mg/L} * V1 = 0.01\text{mg/L} * 19\text{mL}$$

$$V1 = 1.9\text{mL}$$

$$19\text{mL of } 0.01\text{mg/L leachate} = 17.1\text{mL of FSW} + 1.9\text{mL of } 0.1\text{mg/L stock}$$

ZFIX

Going from 18.5% to 4% in a ~200mL final volume!

We want a total of 200 mL volume

$$0.185 * V1 = 0.04 * 200\text{mL}$$

V1 = 43.24mL

To get 200mL of 4% ZFIX: Add 43.24 mL of 18.5% ZFIX to 156.76mL of millipore water

Freeze (& shield)

To preserve embryos for RNA extraction and sequencing we will first submerge them in Zymo DNA/RNA Shield, and then promptly transfer them to a -80C freezer. Here, we won't be dunking them in LN2. The combined DNA/RNA Shield and the quick freeze (due to small size) they will experience in the -80 are together sufficient for preserving RNA.

We will first very gently use a tip-clipped transfer pipette to transfer embryos to our 0.5mL screw-cap cryo-tubes. We will then carefully decant any excess FSW out of the tube, and then add 500uL of Zymo RNA/DNA Shield. Gently invert cryo-tubes a few times to ensure embryos are fully submerged. Once cryo-tubes are filled with their samples, arrange them in wax paper boxes and promptly shift them to -80°C freezer for storage until samples are ready to be processed.

Spawn Night

19:00 - 19:30

- ☐ Move each colony to a numbered/named cambro chamber bin or 5 gallon bucket
- ☐ Ensure each colony has enough headspace so that bundles can rise to the surface
- ☐ Reduce blue holding tank water level to below the lip of the chambers
- ☐ Decant water in chambers such that bouyant bundles don't spill out

19:30 - 20:00

- ☐ Organize the following collection materials at the spawning tubs:
 - ☐ 50mL falcon tubes (1 for ea. colony) filled with 30mL of 0.22micron FSW
 - ☐ Tip-clipped transfer pipettes
 - ☐ Red Headlamps (new batteries)
 - ☐ Printed waterproof spawn collection data sheet
- ☐ Organize the following collection materials at the dry lab bench:
 - ☐ 20mL scintillation vials prepped with treatments
 - ☐ 50mL falcon tube rack (for working with active bundle-bundle crosses)
 - ☐ Tip-clipped transfer pipettes (for moving bundles to scintillation vials)
 - ☐ Printed waterproof bundle-bundle cross metadata sheet (for recording parent colonies of crosses)

20:50 - 21:20

- ☐ Observe for setting (polyp expansion and bundle visibility) via red headlamp
- ☐ Find a spawner!
- ☐ Collect 24 bundles from each spawning colony and transfer them to a pre-labeled 50mL falcon tube with 30mL of FSW
- ☐ Once we have at least 2 colonies worth of bundles, 1 person will move with any filled falcon tubes to the dry lab bench and begin bundle-bundle crosses

- ☐ 1 person will remain and continue to collect bundles from colonies, periodically transferring them to the bundle-bundle crossing station at the dry lab bench

21:20 - 22:00

- ☐ Cross fertilize by taking one bundle each from two distinct colonies and transferring them to the same scintillation vial
- ☐ Record the parent colonies of each cross on the embryonic development metadata sheet
- ☐ Only transfer intact bundles! Once eggs break up and begin to hydrate, the spawn party is over

22:00 - 22:30

- ☐ Observe scintillation vials for bundle breakup and egg hydration
- ☐ Record the times when the first and last vials experience hydration
 - ☐ First Egg Hydration Time _____
 - ☐ Last Egg Hydration Time _____
- ☐ Average the times above , and add 4 hours, 9 hours, and 14 hours
 - ☐ Average Egg Hydration Time _____
 - ☐ 4 hours post-fertilization ¹ _____
 - ☐ 9 hours post-fertilization _____
 - ☐ 14 hours post fertilization _____
- Set alarms for ~40 minutes prior to the times above

- [] (& on July 6th make another batch of Teachate!)

22:30 - 23:00

- ☐ Return each colony to blue tank and remove chambers
- ☐ Return water height to normal
- ☐ Clean up anything left outside

23:00 PM - 01:45

- ☐ 2.5 hr nap time!

02:15 - 03:30

- ☐ ZFIX 4hpf samples (4Z)
 - ☐ Use a transfer pipette to gently move embryos from scintillation vials to microcentrifuge tubes
 - ☐ Decant any excess FSW with transfer pipette
 - ☐ Add ~1mL of 4% ZFIX to each microcentrifuge tube in the fume hood!
 - ☐ Gently invert each capped scintillation vial to mix the ZFIX into the filtered seawater and sample
 - ☐ Record time of fixin' : _____

- ☐ Add 8-12 hours to fixin' time to transfer from ZFIX to 70% ethanol: _____
- ☐ Set fixed vials in 4C fridge
- ☐ Freeze 4hpf samples (4F)
 - ☐ Gently transfer live embryos to cryo-vials
 - ☐ Decant excess FSW using a transfer pipette
 - ☐ Add 500uL of Zymo DNA/RNA Shield
 - ☐ Promptly transfer to -80C freezer
- ☐ Make PVC Leachate (on July 6th!)

03:30 - 06:45

- ☐ 3.25 hr nap time!

07:15 - 08:30

- ☐ ZFIX 9hpf samples (9Z)
 - ☐ Use a transfer pipette to gently move embryos from scintillation vials to microcentrifuge tubes
 - ☐ Decant any excess FSW with transfer pipette
 - ☐ Add ~1mL of 4% ZFIX to each microcentrifuge tube in the fume hood!
 - ☐ Gently invert each capped scintillation vial to mix the ZFIX into the filtered seawater and sample
 - ☐ Record time of fixin' : _____
 - ☐ Add 8-12 hours to fixin' time to transfer from ZFIX to 70% ethanol: _____
 - ☐ Set fixed vials in 4C fridge
- ☐ Freeze 9hpf samples (9F)
 - ☐ Gently transfer live embryos to cryo-vials
 - ☐ Decant excess FSW using a transfer pipette
 - ☐ Add 500uL of Zymo DNA/RNA Shield
 - ☐ Promptly transfer to -80C freezer

08:30 - 11:45

- ☐ 3.25 hr nap time!

12:15 - 13:30

- ☐ ZFIX 14hpf samples (14Z)
 - ☐ Use a transfer pipette to gently move embryos from scintillation vials to microcentrifuge tubes
 - ☐ Decant any excess FSW with transfer pipette
 - ☐ Add ~1mL of 4% ZFIX to each microcentrifuge tube in the fume hood!
 - ☐ Gently invert each capped scintillation vial to mix the ZFIX into the filtered seawater and sample

- ☐ Record time of fixin' : _____
- ☐ Add 8-12 hours to fixin' time to transfer from ZFIX to 70% ethanol: _____
- ☐ Set fixed vials in 4C fridge
- ☐ Freeze 14hpf samples (14F)
 - ☐ Gently transfer live embryos to cryo-vials
 - ☐ Decant excess FSW using a transfer pipette
 - ☐ Add 500uL of Zymo DNA/RNA Shield
 - ☐ Promptly transfer to -80C freezer

Ethanol Transfer Times

4hpf: _____

9hpf: _____

14hpf: _____

::: callout-tip

You can do this!

:::