

Spawn night protocol

For embryonic development of Montipora capitata rice corals

AUTHOR
Sarah Tanja

PUBLISHED
March 26, 2024

Prep

Make Stock Leachate

First make 400mL of each stock solution from the 1000 mg/L prepared leachate

☐ Make 10 mg/L stock:

$$1000 \text{ mg/L} * V1 = 10 \text{ mg/L} * 400\text{mL}$$

$$V1 = 4\text{mL}$$

$$400\text{mL of } 10\text{mg/L stock} = 4\text{mL of } 1000\text{mg/L stock} + 396\text{mL of FSW}$$

☐ Make 1mg/L stock:

$$10\text{mg/L} * V1 = 1\text{mg/L} * 400\text{mL}$$

$$V1 = 40\text{mL}$$

$$400\text{mL of } 1\text{mg/L stock} = 40\text{mL of } 10\text{mg/L stock} + 360\text{mL of FSW}$$

☐ Make 0.1 mg/L stock:

$$1\text{mg/L} * V1 = 0.1\text{mg/L} * 400\text{mL}$$

$$V1 = 40\text{mL}$$

$$400\text{mL of } 0.1\text{mg/L stock} = 40\text{mL of } 1\text{mg/L stock} + 360\text{mL of FSW}$$

2. Then dilute the stock to each treatment vial/jar

For a 20mL scintillation vial for coral embryo experiments with a final target volume of 15mL:

1mg/L HIGH treatment

$$10\text{mg/L} * V1 = 1\text{mg/L} * 19\text{mL}$$

Edit Details

Delete

Spawn night protocol by sarahtanja Last updated 1 minute ago

Comments (–)

Share

Hide Toolbars